




# Calorie Answer CA-Hi

Nutrition analysis by NIR technology

---



Calorie Answer



POWER WAITING CONTROL MEASURING FINISH CANCEL



ON/OFF CONTROL

Jwp

# Agenda

1. 제조사 소개
2. Calorie Answer란?
3. Calorie Answer의 원리 및 측정법별 비교
4. Calorie Answer의 구성 및 사용법
5. 신제품(CA-Hi) 소개
6. Calorie Answer의 활용방안
7. Q&A

# 제조사



Company name : Joy World Pacific Co., Ltd

Location : Aomori Prefecture, Japan

Establishment : October, 1978

Business : Manufacture optical device, Software,  
Hardware development

## History

- 1998 - Start to develop a device to measure Calorie
- 2006 - First model CA-HN in the market
- 2011 - Second model CA-HM in the market

# What is Calorie Answer?



Measured  
parameter

1. Measure standard cell
2. Homogenize and grind sample
3. Put the sample in cell and measure
4. Analysis time is 80 sec.


- kcal
- Protein(g)
- Fat(g)
- Carbohydrate(g)
- Moisture(g)
- Alcohol (g/%)

# Calorie Answer의 특징


- ✓ Use Near Infrared Light
- ✓ Quick and easy measurement operation
- ✓ Possible to measure nutrition of food without knowing those recipes
- ✓ Nondestructive and Non-contact Inspection
- ✓ Patented Technology




# Calorie Answer의 원리



# Calorie Answer의 원리



# 칼로리 측정법 별 비교



		Calorie (kcal)		Carbohydrate (g)		Protein (g)		Fat (g)	
Marinated Spinach	Calculation	30.4		7.9		3.0		0.4	
	A	34.0	A/B	4.7	A/B	2.9	A/B	0.4	A/B
	B	68.1	100%	11.5	145%	4.4	52%	0.5	25%
Fried Chinese Noodle(Salty)	Calculation	170.4		14.8		5.1		5.7	
	A	168.0	A/B	14.9	A/B	4.5	A/B	8.9	A/B
	B	173.8	4%	18.0	21%	6.1	36%	4.4	14%
Spaghetti Napolitans	Calculation	245.2		35.8		7.5		7.2	
	A	145.0	A/B	23.2	A/B	6.1	A/B	3.1	A/B
	B	181.7	25%	29.4	27%	6.8	12%	4.1	32%



# 칼로리 측정법 별 비교

NDL Home Food Search Nutrients List Ground Be

Enter one or more terms

Restaurant, Chinese

[Try a new search](#)

19 foods found Click on a food name to view details

NDB No.	Description
36603	Restaurant, Chinese, beef and vegetables
36617	Restaurant, Chinese, lemon chicken
36620	Restaurant, Chinese, shrimp and vegetables
36626	Restaurant, Chinese, chicken and vegetables
36629	Restaurant, Chinese, orange chicken
36601	Restaurant, Chinese, egg rolls, assorted
36618	Restaurant, Chinese, general tso's chicken
36619	Restaurant, Chinese, kung pao chicken
36621	Restaurant, Chinese, sweet and sour chicken
36622	Restaurant, Chinese, sweet and sour pork
36623	Restaurant, Chinese, chicken chow mein
36602	Restaurant, Chinese, fried rice, without meat


What ingredient is used?



**Fried rice with meat!!**  
**But no data in database.**



How much oil is absorbed?



# 칼로리 측정법 별 비교

	Time	Sample prep	Accuracy	Cost	Repeatability
Laboratory based analysis	X	X	◎	X	??
Database analysis	△	◎	△	◎	??
Calorie Answer	◎	○	○	○	◎


◎ : Excellent

○ : Good

△ : Poor

X : Bad

# Example to distinguish Sucrose and Glucose



Similar

Breakdown

Food ingredients are divided into individual component information

Food  
component

Preinstalled  
Analytical curve

Multiple  
regression  
analysis

# Measurement mode

Many food type,  
component type, form  
type



Choose best  
measurement mode  
depends on food type

## Preinstalled standard mode

Processed foods	Reflection
Snacks	Reflection
Confectionaries	Reflection
Soft drinks	Penetration
Coffee drinks	Penetration

## Other mode(option)

**Potatoes or starches, Sugars, Beans  
Nut, Vegetables, Mushrooms  
Seaweeds, Fishes, Meats, Eggs**


## Additional mode(option)

Cereals	
Cereals SCREENING(also powdered rice)	Reflection
Cereals(under 200kcal)	Reflection
Powdered Cereals, Dried noodles	Reflection


Dairy product	
Yogurt	Reflection
Cheese	Reflection
Powdered milk	Reflection
Yogurt drink	Penetration
Lactic acid bacteria beverage	Penetration
Ice cream, Condensed milk	Penetration
Milk(w/o sugar added)	Penetration
Milk(w/ sugar added)	Penetration
Milk(Whipping cream)	Penetration

# Test result 1


Calorie




Protein




Fat



炭水化物の相関 n=29



水分の相関 n=20



Carb

Moisture

# Test result 2

Butter	Calorie(kcal)	Protein(g)	Fat(g)	Carb(g)	Moisture(g)
Database	745	0.6	81	0.2	16.2
Ave of CA-HM	712	0	78.5	1.6	14.7
Difference Database vs. CA	4%	N.D	3%	88%	10%
CA Measure 1st	(711)	(0.0)	(78.3)	(1.6)	(14.0)
CA Measure 2nd	(716)	(0.0)	(79.0)	(1.3)	(15.0)
CA Measure 3 <sup>rd</sup>	(711)	(0.0)	(78.2)	(1.8)	(15.0)

Beer	Calorie(kcal)	Protein(g)	Fat(g)	Carb(g)	Moisture(g)
Database	40	0.3	0	3.1	92.8
Ave of CA-HM	39	0	0	3.3	92.0
Difference Database vs. CA	2%	N.D	0%	6%	1%
CA Measure 1st	(37)	(0.0)	(0.0)	(2.8)	(93.0)
CA Measure 2nd	(37)	(0.1)	(0.0)	(2.8)	(92.0)
CA Measure 3rd	(43)	(0.0)	(0.0)	(4.2)	(91.0)

# Test result 3

Biscuit	Calorie(kcal)	Protein(g)	Fat(g)	Carb(g)	Moisture(g)
Database	432	7.6	10.0	77.8	2.6
Ave of CA-HM	434	8.3	12.9	71.0	3.0
Difference Database vs. CA	0.5%	8%	23%	10%	14%
CA Measure 1st	(434)	(8.4)	(12.8)	(71.3)	(3.0)
CA Measure 2nd	(434)	(8.4)	(12.9)	(71.1)	(3.0)
CA Measure 3 <sup>rd</sup>	(433)	(8.2)	(13.1)	(70.7)	(3.0)

Bread	Calorie(kcal)	Protein(g)	Fat(g)	Carb(g)	Moisture(g)
Database	264	9.3	4.4	46.7	38.0
Ave of CA-HM	286	11.8	8.8	39.9	38.7
Difference Database vs. CA	8%	21%	50%	17%	2%
CA Measure 1st	(286)	(12.0)	(8.8)	(39.7)	(38.0)
CA Measure 2nd	(286)	(11.7)	(8.9)	(39.9)	(39.0)
CA Measure 3 <sup>rd</sup>	(286)	(11.6)	(8.8)	(40.0)	(39.0)

# System configuration



PC is not included in the system.





# How to sample preparation





Homogenize sample

Reflection measurement



Penetration measurement



# Sample analysis

Reference cell measurement



Sample measurement



Reflection measurement

Measure  
without  
cell



Penetration measurement

Only  
5 min for  
analysis!!

# Output data

JWP CA-HM Client 2.5.5

File Edit Help

Icons: [Print] [Print with Selection] [Zoom In] [Zoom Out] [Close] [Close All]

Sample Name	Measure...	Tempera...	Weight (g)	Energy (kcal)	Protein (g)	Fat (g)	Carbohyd...	S
M11検証 オリーブオ...	油脂類 (...)	16.4	100	900	0.0	100.0	0.0	
M11検証 オリーブオ...	油脂類 (...)	16.6	100	900	0.0	99.9	0.1	
M11検証 オリーブオ...	油脂類 (...)	16.6	100	900	0.0	100.0	0.0	
M11検証 オリーブオ...	油脂類 (...)	16.6	100	900	0.0	100.0	0.0	
M11検証 オリーブオ...	油脂類 (...)	16.3	100	900	0.0	100.0	0.0	
M11検証 ごま油 (環...	油脂類 (...)	16.4	100	900	0.0	100.0	0.0	
M11検証 ごま油 (環...	油脂類 (...)	15.8	100	900	0.0	100.0	0.0	

# Improvement CA-Hi

## 1. Reference cell measurement

- Reference cell은 최초 1회만 측정→총 측정시간 감소  
(최초 측정 약 80초, 후속 측정 약 80초 소요)

## 2. Analysis of Liquid sample

- “transmission reflection method(투과반사법)”을 사용  
→액체가 내부로 새는 것을 방지

## 3. Analysis mode


- 기존장비: 6가지 측정모드 선택, 설치  
→calibration curve creation software 설치가능  
multiple measurement modes(측정모드에 제한이 없음)

# New Software

3. 波形表示  
(파형데이터표시)

5. 検量線作成  
(분석곡선작성)

7. 測定値  
(측정치)



# What the system cannot do




Zero-Calorie



Black Color Foods

Count calorie of unavailable carb due to molecule that consists of sugar.


Black color food absorb light well results in poor reflection light.



# Low digestible carb

## Food containing unavailable carb, such as erythritol

**Unavailable carb**  
**Is low digestion by human**  
**such as sugar alcohols or**  
**Artificial sweeteners**



Erythritol : 0kcal/g  
Agar : 2kcal/g  
Maltitol : 3kcal/g  
Sorbitol : 3kcal/g  
Xylitol




**Use Erythritol (14g)**




Calorie : 352 kcal  
Protein : 6.2g  
Fat : 24.1g  
Carb : 27.5g  
Moisture : 41%

14g of Erythritol is  
contained in 27.5g of Carb

  
**Manual  
calculation**

**Available calorie for human is**  
**296 kcal**  
**(352 – 14 X 4)**

# Use for



- ✓ Display of calorie and nutrition value
- ✓ New product, new recipe developments
- ✓ Quality Control
- ✓ Nutrition management by checking leftover foods



# Use for new products development

## Cooking method, ingredients for low calorie



- Ingredients?
- Cooking method?



- What kind of oil?
- How much oil absorbed?

# Use for Calorie management



- Precise calorie management for dieter, athlete
- Unique concept restaurant, training gym



Spaghetti Napolitano  
12kcal



Hamburg steak  
15kcal




Low calorie restaurant 『etsu』

In general, calorie restricted diet is considered “unsavory”. This restaurant serves low calorie, savory dish. They also sell unique food stuff for making low calorie food in home.



# Calorie answer User



# Best for

## Purpose

## Object

## Effectiveness

- Display of calorie
- Insistence of company stance
- Increase efficiency of dietitian
- Substitution of physicochemical analysis

- Restaurant chain
- Daily foods shops

- Increase of consumer sufficient
- Increase of marquee
- Cost down

- New product development
- Development of cooking devices
- Manufacturing process control
- Control of raw materials

- Restaurant Chain
- Food industries
- Electronic industries
- Cooking devices maker

# Contact info

## 씨엔티교역

(T)031-308-9777, (F)031-308-9779

Email: [tastesenser@naver.com](mailto:tastesenser@naver.com),

[foodscience7@naver.com](mailto:foodscience7@naver.com),

[cnt@cntglobal.co.kr](mailto:cnt@cntglobal.co.kr)